

Accurate Speed Control

FMI RH-ASC / BM BM-FMI QP-pumps







Oina VV AB SE-139 90 Stäksholmen Värmdö (Stockholm) Sweden Tel: +46 8 410 163 70 info@oinapumps.com www.oinapumps.com

USERS MANUAL	3
WHAT IS IN THE BOX	3
HOW TO SET UP THE ACCURATE SPEED CONTROL BM / FMI SYSTEM	
FUNCTIONALITY IN MODE	4
How to set MODE functions	4
REMOTE CONTROL FUNCTIONS VIA 15 PIN D-SUB	5
External Speed Control	5
REAR PANEL CONNECTIONS	6

User Guide

What is in the box

- 1 ea Pump Unit with inbuilt electronic for DC motor control and timer and FMI pump head mounted on the front panel.
- (24VDC, 30W)
- 1 ea Power Supply (Input 100-240VAC, 50/60Hz, 800mA. Output 24VDC, 35VA)
- 1 ea User Guide
- Customised remote cables, foot pedals or other accessories may added on request

WARNING!

No other parts will be used for this system except parts delivered from Oina VV AB! Oina VV AB is not responsible for any injury or damage caused by misconduct!

How to set up the ACCURATE Speed Control BM / FMI system

- Connect the enclosed Power Supply to the Pump Unit, and then connect the Power Supply to the main power, 100VAC or 230VAC, (100-240VAC, 50/60Hz, 800mA)
- Connect remote control cable or foot pedal or other accessories if requested
- First time:
 - <u>Factory default setting:</u>
 - Red LED with **STOP** is switched on.
 - Green LED with "INT-SPEED" switched on.
 - 50 % Speed is chosen.
- Now you are ready to run the pump as a continuous metering pump.
 - Press **FWD** button to run in forward direction or press **REV** button to run in reverse direction.
 - \circ Adjust speed from 5-100% in steppes of 1% via + and press buttons.
 - \circ For momentary flush the pump, press **MAX** button.
- The Pump Unit will remember the last setting when it will be reconnected next time to main power or after power failure.

Functionality in MODE

- **TIMER:** 1-999 seconds in steppes of 1 second
- **COUNTER:** Not available function on this Pump Unit!
- INTERNAL SPEED CONTROL
- **EXTERNAL SPEED CONTROL:** For communication with process instrumentation via0-20mA or 4-20mA
- **I/O CONFIGURATION:** Set up for RS232. This function is not available on this Pump Unit!
- **RATIO/DIR:** Factory set up configurations for different motors, gear ratios and for main direction. Users are not able to make changes!

How to set MODE functions

• Timer-Dispensing

Press **MODE** button multiple times until green LED lights up with TIMER.

(First time you press **MODE** the red LED will be switched off!)

Set dispensing time via + and - buttons from 1-999 seconds.

Press **STOP** to confirm your choice of dispensing time.

Red LED with **STOP**, and the green LED's with "TIMER" and with "INT-SPEED" will be switched on to confirm your choice.

Adjust speed from 5-100% in steppes of 1% via + and - press buttons.

Press **FWD** button for each dispensing in forward direction or press **REV** button for each dispensing in reverse direction.

• To leave Dispensing and go back to Metering function

Press MODE button multiple times until green LED lights up with "INT-SPEED" (or "EXT-SPEED" if external speed control is chosen). Set dispensing time to 000 Press STOP button to confirm. Red LED with STOP button, and green LED's with "INT-SPEED" (or "EXT-SPEED" if external speed control is chosen) will switch on to confirm your choice. Adjust speed from 5-100% in steppes of 1% via + and - press buttons. Start metering via FWD or REV press buttons.

Remote Control Functions via 15 pin D-sub

External Speed Control

Press MODE button multiple times until green LED lights up with "EXT-SPEED".
Scroll in-between - - -, 0.20 and 4.20 via + and - press buttons:
- - means "INT-SPEED" (Factory set)
0.20 means external speed control via 0-20mA signal
4.20 means external speed control via 4-20mA signal
When right choice shows, press MODE button to confirm.

External Start

Run FWD or run Rev Start dispensing FWD or start dispensing REV via customised cable, switch or foot pedal!

Counter function and RS232 via 9 pin D-sub, for this pump is not in function

For external speed control, see users manual attached to customised cable!

For security reasons only cables or foot pedals delivered by Oina VV AB are allowed to be connected to our products! For all kind of cables or foot pedals connected to the Pump Unit for external speed control, start or dispense forward or reverse: read users manual attached to customised foot pedal or customised remote cable

Accurate Speed Control accommodates all FMI pump head sizes.

See separate users manual for all info about the FMI pump heads!

Tel: +46 8 410 163 70 **www.oinapumps.com**

Rear panel connections

