

# Automatic Liquid Handling - Stroke 60 mm

**AD**

CGA\* EPOS / EASY / ERIS, Roche\* Midias



## TECHNICAL INFORMATION

Glass Barrel: Borosilicate 3.3

### PTFE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Part No.
50 µL	H	60 mm	PTFE	M8x0.75	2605613
100 µL	H	60 mm	PTFE	M8x0.75	2605713
250 µL	H	60 mm	PTFE	M8x0.75	2605813
500 µL	H	60 mm	PTFE	M8x0.75	2605913
1.0 mL	H	60 mm	PTFE	M8x0.75	2605014
2.5 mL	H	60 mm	PTFE	M8x0.75	2605034
5.0 mL	H	60 mm	PTFE	M8x0.75	2605054
10.0 mL	H	60 mm	PTFE	M8x0.75	2605074
25.0 mL	H	60 mm	PTFE	M8x0.75	2605094

### PTFE-Seal with Silicone

Volume	Series	Stroke	Plunger Tip	Front Fitting	Part No.
50 µL	H	60 mm	Silicone	M8x0.75	2605617
100 µL	H	60 mm	Silicone	M8x0.75	2605717
250 µL	H	60 mm	Silicone	M8x0.75	2605817
500 µL	H	60 mm	Silicone	M8x0.75	2605917
1.0 mL	H	60 mm	Silicone	M8x0.75	2605017
2.5 mL	H	60 mm	Silicone	M8x0.75	2605037
5.0 mL	H	60 mm	Silicone	M8x0.75	2605057
10.0 mL	H	60 mm	Silicone	M8x0.75	2605077

### HDPE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Part No.
50 µL	H	60 mm	HDPE	M8x0.75	2605618
100 µL	H	60 mm	HDPE	M8x0.75	2605718
250 µL	H	60 mm	HDPE	M8x0.75	2605818
500 µL	H	60 mm	HDPE	M8x0.75	2605918
1.0 mL	H	60 mm	HDPE	M8x0.75	2605018
2.5 mL	H	60 mm	HDPE	M8x0.75	2605038
5.0 mL	H	60 mm	HDPE	M8x0.75	2605058
10.0 mL	H	60 mm	HDPE	M8x0.75	2605078

\*All other trademarks, service marks, registered trademarks, and registered service marks are the property of their respective owners. Gardner Denver and ILS are not affiliated with CGA and Roche in any way.