

Automatic Liquid Handling - Stroke 30 mm

H-PD



TECHNICAL INFORMATION

Glass Barrel: Borosilicate 3.3

PTFE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
12.5 µL	H	30 mm	PTFE	Luer	6-32 UNC	2658415
25 µL	H	30 mm	PTFE	Luer	6-32 UNC	2658515
50 µL	H	30 mm	PTFE	Luer	6-32 UNC	2658615
100 µL	H	30 mm	PTFE	Luer	6-32 UNC	2658715
125 µL	H	30 mm	PTFE	Luer	6-32 UNC	2658714
250 µL	H	30 mm	PTFE	Luer	6-32 UNC	2658815
500 µL	H	30 mm	PTFE	Luer	6-32 UNC	2658915
1.0 mL	H	30 mm	PTFE	Luer	6-32 UNC	2658014
1.25 mL	H	30 mm	PTFE	Luer	6-32 UNC	2658015
2.5 mL	H	30 mm	PTFE	Luer	6-32 UNC	2658035
5.0 mL	H	30 mm	PTFE	Luer	6-32 UNC	2658055
12.5 mL	H	30 mm	PTFE	Luer	6-32 UNC	2658074

PTFE-Seal with Silicone

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
12.5 µL	H	30 mm	Silicone	Luer	6-32 UNC	2666415
25 µL	H	30 mm	Silicone	Luer	6-32 UNC	2666515
50 µL	H	30 mm	Silicone	Luer	6-32 UNC	2666615
100 µL	H	30 mm	Silicone	Luer	6-32 UNC	2666715
125 µL	H	30 mm	Silicone	Luer	6-32 UNC	2666714
250 µL	H	30 mm	Silicone	Luer	6-32 UNC	2666815
500 µL	H	30 mm	Silicone	Luer	6-32 UNC	2666915
1.0 mL	H	30 mm	Silicone	Luer	6-32 UNC	2666014
1.25 mL	H	30 mm	Silicone	Luer	6-32 UNC	2666015
2.5 mL	H	30 mm	Silicone	Luer	6-32 UNC	2666035
5.0 mL	H	30 mm	Silicone	Luer	6-32 UNC	2666055
12.5 mL	H	30 mm	Silicone	Luer	6-32 UNC	2666074

UHMWPE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
25 µL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665515
50 µL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665615
100 µL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665715
125 µL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665714
250 µL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665815
500 µL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665915
1.0 mL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665014
1.25 mL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665015
2.5 mL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665035
5.0 mL	H	30 mm	UHMWPE	Luer	6-32 UNC	2665055