


# Automatic Liquid Handling - Stroke 60 mm

## ILS Syringe compatible with ABI PRISM 3100 Genetic Analyzer



### TECHNICAL INFORMATION

Glass Barrel: Borosilicate 3.3

#### Fitting w/o insert

Volume	Plunger Tip	Front Fitting	Part No.
250 µL	PTFE	1/4"-32 UNEF w/o insert	2618811
500 µL	PTFE	1/4"-32 UNEF w/o insert	2618911
1.0 mL	PTFE	1/4"-32 UNEF w/o insert	2618011
5.0 mL	PTFE	1/4"-32 UNEF w/o insert	2618051

#### Fitting with insert

Volume	Plunger Tip	Front Fitting	Part No.
250 µL	PTFE	1/4"-32 UNEF with insert	2618812
500 µL	PTFE	1/4"-32 UNEF with insert	2618912
1.0 mL	PTFE	1/4"-32 UNEF with insert	2618012
5.0 mL	PTFE	1/4"-32 UNEF with insert	2618052
10.0 mL	PTFE	1/4"-32 UNEF with insert	2618081