

ILS Syringe compatible with DAD (Diluter ACCUDIL)



TECHNICAL INFORMATION

Side Connection: M6x1

Glass Barrel: Borosilicate 3.3

PTFE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
50 µL	H	60 mm	PTFE	M8x0.75	6-32 UNC	2607618
100 µL	H	60 mm	PTFE	M8x0.75	6-32 UNC	2607718
250 µL	H	60 mm	PTFE	M8x0.75	6-32 UNC	2607818
500 µL	H	60 mm	PTFE	M8x0.75	6-32 UNC	2607918
1.0 mL	H	60 mm	PTFE	M8x0.75	6-32 UNC	2607018
2.5 mL	H	60 mm	PTFE	M8x0.75	6-32 UNC	2607038

HDPE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
50 µL	H	60 mm	HDPE	M8x0.75	6-32 UNC	2626618
100 µL	H	60 mm	HDPE	M8x0.75	6-32 UNC	2626718
250 µL	H	60 mm	HDPE	M8x0.75	6-32 UNC	2626818
500 µL	H	60 mm	HDPE	M8x0.75	6-32 UNC	2626918
1.0 mL	H	60 mm	HDPE	M8x0.75	6-32 UNC	2626018
2.5 mL	H	60 mm	HDPE	M8x0.75	6-32 UNC	2626038

PTFE-Seal with Silicone

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
50 µL	H	60 mm	Silicone	M8x0.75	6-32 UNC	2626617
100 µL	H	60 mm	Silicone	M8x0.75	6-32 UNC	2626717
250 µL	H	60 mm	Silicone	M8x0.75	6-32 UNC	2626817
500 µL	H	60 mm	Silicone	M8x0.75	6-32 UNC	2626917
1.0 mL	H	60 mm	Silicone	M8x0.75	6-32 UNC	2626017
2.5 mL	H	60 mm	Silicone	M8x0.75	6-32 UNC	2626037

