


## Automatic Liquid Handling - Stroke 60 mm

### ILS Syringe compatible with DIONEX AS



#### TECHNICAL INFORMATION

Glass Barrel: Borosilicate 3.3

#### PTFE-Seal

Volume	Stroke	Plunger Tip	Front Fitting	Part No.
25 µL	60 mm	PTFE	1/4"-28 UNF	2608525
50 µL	60 mm	PTFE	1/4"-28 UNF	2608625
100 µL	60 mm	PTFE	1/4"-28 UNF	2608725
250 µL	60 mm	PTFE	1/4"-28 UNF	2608825
500 µL	60 mm	PTFE	1/4"-28 UNF	2608925
1.0 mL	60 mm	PTFE	1/4"-28 UNF	2608025

#### HDPE-Seal

Volume	Stroke	Plunger Tip	Front Fitting	Part No.
25 µL	60 mm	HDPE	1/4"-28 UNF	2608527
50 µL	60 mm	HDPE	1/4"-28 UNF	2608627
100 µL	60 mm	HDPE	1/4"-28 UNF	2608727
250 µL	60 mm	HDPE	1/4"-28 UNF	2608827
500 µL	60 mm	HDPE	1/4"-28 UNF	2608927
1.0 mL	60 mm	HDPE	1/4"-28 UNF	2608027