


# Automatic Liquid Handling - Stroke 60 mm

## RQ

Zinsser Analytic Lissy / Sophas / Speedy / Calli / Moss, Rosys\* Plato & Lucy, Qiagen\* Biorobot 9600, Digene\* Rapid Capture System, Unitech Chemwell\* 2910



## TECHNICAL INFORMATION

Glass Barrel: Borosilicate 3.3

### PTFE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
25 µL	H	60 mm	PTFE	Luer	Screw	2607512
50 µL	H	60 mm	PTFE	Luer	Screw	2607612
100 µL	H	60 mm	PTFE	Luer	Screw	2607712
250 µL	H	60 mm	PTFE	Luer	Screw	2607812
500 µL	H	60 mm	PTFE	Luer	Screw	2607912
1.0 mL	H	60 mm	PTFE	Luer	Screw	2607014
2.5 mL	H	60 mm	PTFE	Luer	Screw	2607034
5.0 mL	H	60 mm	PTFE	Luer	Screw	2607054
10.0 mL	H	60 mm	PTFE	Luer	Screw	2607074

### PTFE-Seal with Silicone

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
25 µL	H	60 mm	Silicone	Luer	Screw	2644517
50 µL	H	60 mm	Silicone	Luer	Screw	2644617
100 µL	H	60 mm	Silicone	Luer	Screw	2644717
250 µL	H	60 mm	Silicone	Luer	Screw	2644817
500 µL	H	60 mm	Silicone	Luer	Screw	2644917
1.0 mL	H	60 mm	Silicone	Luer	Screw	2644017
2.5 mL	H	60 mm	Silicone	Luer	Screw	2644037
5.0 mL	H	60 mm	Silicone	Luer	Screw	2644057
10.0 mL	H	60 mm	Silicone	Luer	Screw	2644077

### HDPE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Plunger Top	Part No.
50 µL	H	60 mm	HDPE	Luer	Screw	2644618
100 µL	H	60 mm	HDPE	Luer	Screw	2644718
250 µL	H	60 mm	HDPE	Luer	Screw	2644818
500 µL	H	60 mm	HDPE	Luer	Screw	2644918
1.0 mL	H	60 mm	HDPE	Luer	Screw	2644018
2.5 mL	H	60 mm	HDPE	Luer	Screw	2644038
5.0 mL	H	60 mm	HDPE	Luer	Screw	2644058
10.0 mL	H	60 mm	HDPE	Luer	Screw	2644078