

# Automatic Liquid Handling - Stroke 30 mm

**XE**



## TECHNICAL INFORMATION

Glass Barrel: Borosilicate 3.3

### PTFE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Part No.
25 µL	H	30 mm	PTFE	1/4"-28 UNF	2607521
50 µL	H	30 mm	PTFE	1/4"-28 UNF	2607621
100 µL	H	30 mm	PTFE	1/4"-28 UNF	2607721
250 µL	H	30 mm	PTFE	1/4"-28 UNF	2607821
500 µL	H	30 mm	PTFE	1/4"-28 UNF	2607921
1.0 mL	H	30 mm	PTFE	1/4"-28 UNF	2607021
1.5 mL	H	30 mm	PTFE	1/4"-28 UNF	2607027
2.5 mL	H	30 mm	PTFE	1/4"-28 UNF	2607041
5.0 mL	H	30 mm	PTFE	1/4"-28 UNF	2607061

### PTFE-Seal with Silicone

Volume	Series	Stroke	Plunger Tip	Front Fitting	Part No.
500 µL	H	30 mm	Silicone	1/4"-28 UNF	2646927
1.0 mL	H	30 mm	Silicone	1/4"-28 UNF	2646027
2.5 mL	H	30 mm	Silicone	1/4"-28 UNF	2646047
5.0 mL	H	30 mm	Silicone	1/4"-28 UNF	2646067

### HDPE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Part No.
50 µL	H	30 mm	HDPE	1/4"-28 UNF	2646628
100 µL	H	30 mm	HDPE	1/4"-28 UNF	2646728
500 µL	H	30 mm	HDPE	1/4"-28 UNF	2646928
1.0 mL	H	30 mm	HDPE	1/4"-28 UNF	2646028
2.5 mL	H	30 mm	HDPE	1/4"-28 UNF	2646048
5.0 mL	H	30 mm	HDPE	1/4"-28 UNF	2646068