

Automatic Liquid Handling - Stroke 30 mm

XP-PD



TECHNICAL INFORMATION

Plunger Top: 6-32UNC

Glass Barrel: Borosilicate 3.3

PTFE-Seal

Volume	Series	Stroke	Plunger Tip	Front Fitting	Part No.
12.5 µL	H	30 mm	PTFE	1/4"-28 UNF	2676401
25 µL	H	30 mm	PTFE	1/4"-28 UNF	2676501
50 µL	H	30 mm	PTFE	1/4"-28 UNF	2676601
100 µL	H	30 mm	PTFE	1/4"-28 UNF	2676701
125 µL	H	30 mm	PTFE	1/4"-28 UNF	2676702
250 µL	H	30 mm	PTFE	1/4"-28 UNF	2676801
500 µL	H	30 mm	PTFE	1/4"-28 UNF	2676901
1.0 mL	H	30 mm	PTFE	1/4"-28 UNF	2676011
1.25 mL	H	30 mm	PTFE	1/4"-28 UNF	2676012
2.5 mL	H	30 mm	PTFE	1/4"-28 UNF	2676031
5.0 mL	H	30 mm	PTFE	1/4"-28 UNF	2676051
12.5 mL	H	30 mm	PTFE	1/4"-28 UNF	2676072