

CYCLOSAVER SAV5 DELIVERY

RPM	PSIG	0.0	7.3	14.5	22	29	36
1000	CFM	254	241	230	218	207	197
1200	CFM	303	292	281	271	262	250
1400	CFM	360	349	337	328	317	304
1600	CFM	421	406	391	379	368	357

CYCLOSAVER SAV8 DELIVERY

RPM	PSIG	0.0	7.3	14.5	22	29	32
1000	CFM	347	334	323	312	299	294
1200	CFM	424	412	400	388	377	365
1400	CFM	496	486	477	468	456	452
1600	CFM	577	567	556	544	534	530
1800	CFM	655	642	631	619	607	603

Actual performance may vary based on system design, environment, and maintenance.