

Pump Model CCE10

The Solution to Your Application

Gardner Denver model MTO II is changing the way the world uses oilfree compressed air. By developing a compressor that can run 24 hours a day, Gardner Denver has eliminated the need for over-sizing and back-up units that add unnecessary expense. The craftsmanship of the MTO II is nothing short of the best in the industry. Whether it's healthcare, instrumentation, or laboratory applications, the need for pure, environmentally clean compressed air is rapidly expanding. Gardner Denver is the proven choice for oilfree air.

Standard Features

All MTO II compressors are engineered as 100% duty cycle, V-belt-driven packaged units with:

- Long-lasting cast iron, hard chrome-plated cylinders with PTFE guide rings
- Reed valve design
- Strong ductile iron crankshafts
- Large aluminum cooling fan and low compressor speeds
- Inlet air filter silencers
- Enclosed belt guards
- Electric motor
- Discharge line unloader solenoid valves

- Pressure relief valves
- ASME coded receivers
- NFPA Level 1 and 3 compliant systems
- Grade "D" breathing air packages are available
- 1-Year package warranty
- 1-Year pump warranty



BARE PUMP

MODEL	APPROX. SHIP WEIGHT LBS.	BORE DIA. INCHES	STROKE INCHES	DISPLACEMENT CFM MAX @ 100 PSI	COOLING AIR FLOW MAX CFM @ 100 PSI	HEAT REJECTION MAX BTU/HR	MAX PUMP RPM
CCE10	47	2.95	2.17	4.9	450	2100	575

SIMPLEX BASE MOUNTED UNITS

MOTOR	MTO II MATRIX	DIMENSIONS L × W × H	APPROX. SHIP WEIGHT	CFM DEL'Y		
HP	CHBMOA	INCHES	LBS.	50 PSIG	100 PSIG	
0.75	B.75MTOII	36 × 22 × 25	132	2.4	1.7	
1	B1MTOII	36 × 22 × 25	133	3.8	2.7	

SIMPLEX TANK MOUNTED UNITS

MOTOR	MTO II MATRIX	TANK CAPACITY GAL.	DIMENSIONS L × W × H INCHES	APPROX. SHIP WEIGHT LBS.	CFM DEL'Y	
HP	CHSMOA				50 PSIG	100 PSIG
0.75	H.75MTOII-3	30	38 × 25.5 × 37.5	202	2.4	1.7
1	H1MTOII-3	30	38 × 25.5 × 37.5	203	3.8	2.7
1	H1MTOII-6	60	48 × 25.5 × 43	322	3.8	2.7

DUPLEX TANK MOUNTED UNITS

MOTOR HP	MTO II MATRIX CHDMOA	TANK CAPACITY GAL.	DIMENSIONS L × W × H INCHES	APPROX. SHIP WEIGHT LBS.	CFM DEL'Y	
					50 PSIG	100 PSIG
0.75	H.75MTOIID-6	60	70 × 29 × 43	459	4.8	3.4
1	H1MTOIID-6	60	70 × 29 × 43	460	7.6	5.4
1	H1MTOIID-8	80	70 × 29 × 43	520	7.6	5.4



Gardner Denver, Inc.

1800 Gardner Expressway Quincy, IL 62305 866-440-6241

www.gardnerdenver.com/gdproducts









