## NASH 2BE1-7 LIQUID RING COMPRESSOR SERIES



## RELIABLE PERFORMANCE. HIGHER DISCHARGE PRESSURE.

Expanding the pressure performance of the proven reliable, and versatile NASH 2BE1 liquid ring compressor.




The proven reliable, and versatile NASH 2BE1 liquid ring compressor series now reaches higher pressures.

- > Proven reliable NASH 2BE1 performance
- > Higher capacity capability for discharge pressures up to 35 psig (3.5 bar abs.)
- > Efficient performance
- > Large differential pressure capabilities
- Optional mechanical seals (single or double acting with self - or external flushing)
- > Optional O-Ring sealing possible
- > Optional inlet/outlet manifold
- Available in stainless steel, cast iron and other materials
- > High hydro test pressure capability (up to 10 bar with O-Ring sealing)
- > ATEX Category 1 and Category 2
- > NASH Certified™ 2 Year Warranty

## **PERFORMANCE**

CAPACITY	47 - 930 scfm	80 - 1580 Nm³/h
DISCHARGE PRESSURE	35 psig	3.5 bar abs.
INDUSTRIES	Chemical, Petrochemical, General Industrial	
APPLICATIONS	Reaction process, VCM Recovery, Explosive Gas Compression, Filtration, Solvent Recovery, Garoline Vapor Recovery	





by Gardner Denver

GARDNER DENVER NASH, LLC

tel: +1 724 239 1500

nash@gardnerdenver.com

www.GDNash.com

THE ORIGINAL 2BE1 COMPRESSOR SERIES.
AVAILABLE ONLY FROM NASH. SMARTER WORKS HARDER.